

Docman RMS

technical specification

Docman RMS technical specification

Version: 1

Table of contents

Document Information	3
History	3
Copyright.....	3
Exclusion of liability	3
Ownership	3
Introduction	4
Docman RMS database server	5
Minimum specification (50 Clients)	5
Recommended specification (100 Clients).....	5
Recommended specification (250 Clients or more)	5
Docman RMS client machine	6
Minimum specification.....	6
Recommended specification	6

Document Information

History

Author	Stage	Description	Version ¹
Benjamin Rowe	Initial Draft	Creation of first draft.	0.2
Benjamin Rowe	Release	Release document	1.0

Copyright

© The copyright in this document is owned by PCTI Solutions Limited, including (without limitation) in all alpha-numeric text and code, images, schematics and illustrations. No part of this document may be reproduced, stored in a retrieval system, or transmitted in any form or by any means without the express written consent of PCTI Solutions Limited.

Exclusion of liability

The information in this document is provided for information only and is subject to change without notice. While we make strenuous efforts to ensure that the content of this document is accurate and up to date when published, PCTI Solutions Limited gives no warranty on its accuracy and assumes no responsibility or liability for any errors or inaccuracies whatsoever (to the extent such exclusion is permitted by applicable law).

Ownership

For more information on the content or status of this document, please contact the person identified below:

Benjamin Rowe
Technical Services Manager

E-mail: benrowe@docman.com
Tel: 01977 66 44 88

Docman
Pioneer Court
Pioneer Way
Whitwood, Castleford, West Yorkshire, WF10 5QU

Tel: 01977 66 44 88

Introduction

This document outlines the requirements for the PCTI product 'DOCMAN RMS'. It shows the minimum and recommended specifications for hardware and operating system software for a PC based machine.

The Document will outline two areas of specification:

1. The Docman RMS Database Server
2. The Docman client machines

Docman RMS database server

In order to run Docman RMS, a server on which to run the Docman RMS database is required.

Minimum specification (50 Clients)

- Processor speed of 2Ghz (Or 1 Core VCPU)
- 100GB free hard disk space
- 2GB Ram
- Network Card and a network connection
- CD Drive/USB Port
- Microsoft Windows 2008 R2 Server (SP2)
- Internet Explorer 10

Recommended specification (100 Clients)

- Processor Intel Xeon 3.0 GHZ or Equivalent or higher (Or 2 Core VCPU)
- 250GB free hard disk space
- 4GB Ram
- Network Card and a network connection
- CD Drive/USB Port
- Microsoft Windows 2008 R2 Server (SP2)
- Internet Explorer 11

Recommended specification (250 Clients or more)

- Processor Intel Xeon 3.0 GHZ or Equivalent or higher (Or 4 Core VCPU)
- 500GB free hard disk space or more.
- 8GB Ram or more.
- Network Card and a network connection
- CD Drive/USB Port
- Microsoft Windows 2008 R2 Server (SP2)
- Internet Explorer 11

Docman RMS client machine

All workstations requiring access to the Docman RMS database will require a Docman client license. Workstations accessing Docman via a thin-client 'remote desktop' session also require a Docman client license

Each workstation must have network access to the Docman RMS Server machine.

Minimum specification

- Processor speed of 2Ghz
- 10GB free hard disk space
- 1GB Ram
- Network Card and a network connection within the surgery
- CD Drive/USB Port
- Microsoft Windows 7 (SP1)
- Internet Explorer 9

Recommended specification

- Processor Intel Xeon 3.0 GHZ
- 100GB free hard disk space
- 2GB Ram
- Network Card and a network connection within the surgery
- CD Drive/USB Port
- Microsoft Windows 7 (SP1)
- Internet Explorer 11